

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

6-Chloropyridine-2-carboxylic acid - 97%

Product Number:

680044

CAS Number:

4684-94-0

MDL:


MFCD00155390

Formula:

C₆H₄ClNO₂

Formula Weight:

157.55 g/mol

**TEST****Specification**

Appearance (Color)

Conforms to Requirements

White to Grey to Tan

Appearance (Form)

Powder or Crystals or Chunk(s)

Infrared spectrum

Conforms to Structure

Purity (Titration by NaOH)

96.5 - 103.5 %

GC (area %) (VPCT)

≥ 96.5 %

Specification: PRD.0.ZQ5.10000064911

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.