

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name: Adrenosterone - 98%

Product Number: 284998 CAS Number: 382-45-6

MDL: MFCD00003606 Formula: C19H24O3 Formula Weight: 300.39 g/mol H₃C O H H H

TEST Specification

Appearance (Color)

Appearance (Form)

Infrared spectrum Carbon

Optical Rotation

C = 1%, CHCl3

Purity (TLC)

Solubility (Turbidity)

50mg/ml, CHCl3

Solubility (Color)

White

Powder or Solid or Crystals Conforms to Structure

74.1 - 77.9 % 290.0 - 310.0 deg

≥ 98 %

Clear

Colorless to Very Faint Yellow

Specification: PRD.2.ZQ5.10000034131

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.