


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
4-(Trifluoromethoxy)phenylhydrazine hydrochloride - 98%

Product Number: **366951**
CAS Number: 133115-72-7
MDL: MFCD00053033
Formula: C₇H₇F₃N₂O · HCl
Formula Weight: 228.60 g/mol



TEST

Specification

Appearance (Color)
Appearance (Form)
Infrared spectrum
Assay

White to Beige to Brown
Powder or Crystals or Flakes
Conforms to Structure
≥ 97.5 %

Specification: PRD.2.ZQ5.10000054561

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.