

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com

Email USA: techserv@sial.com

Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Product Number: 805386



Formula: C₄₀H₄₅CoF₁₂N₁₁O₈S₄
 Formula Weight: 1,223.02 g/mol

TEST

Specification

Appearance (Color)	Faint Orange to Very Dark Orange
Appearance (Form)	Powder or Crystals
Infrared Spectrum	Conforms to Structure
Carbon Content	38.2 - 40.3 %
Nitrogen Content	12.2 - 13.0 %
Product of Supplier	Confirmed
Greatcell Solar®	
Registered Trademark	Note
Greatcell Solar is a registered trademark of Greatcell Solar	
Recommended Retest Period	-----
2 Years	

Specification: PRD.0.ZQ5.10000097666

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.