


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-Naphthoic acid - 98%

Product Number: 180246
CAS Number: 93-09-4
MDL: MFCD00004101
Formula: C₁₁H₈O₂
Formula Weight: 172.18 g/mol



TEST

Specification

Appearance (Color)

Conforms to Requirements

White to Off-White to Beige

Appearance (Form)

Powder

Infrared spectrum

Conforms to Structure

Purity (Titration by NaOH)

97.5 - 102.5 %

Purity (HPLC)

≥ 97.5 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.