


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
10-Camphorsulfonyl chloride - 98%

Product Number: **301957**  
CAS Number: 4552-50-5  
MDL: MFCD00064155  
Formula: C<sub>10</sub>H<sub>15</sub>ClO<sub>3</sub>S  
Formula Weight: 250.74 g/mol



### TEST

### Specification

Appearance (Color)

Conforms to Requirements

Off-White to Tan or Pale Yellow

Appearance (Form)

Powder

Infrared spectrum

Conforms to Structure

Titration by AgNO<sub>3</sub>

97.5 - 102.5 %

After Oxygen Combustion

Specification: PRD.0.ZQ5.10000027313

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.