

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.com

Email USA: techserv@sial.com


Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

4-Bromobenzenediazonium tetrafluoroborate - 96%

Product Number: 280895
CAS Number: 673-40-5
MDL: MFCD00011894
Formula: C₆H₄BBrF₄N₂
Formula Weight: 270.82 g/mol
Storage Temperature: 2 - 8 °C



TEST

Specification

Appearance (Color)	White to Tan
Appearance (Form)	Conforms to Requirements
Powder, Crystals, Crystalline Powder and/or Chunks	
Infrared spectrum	Conforms to Structure
Carbon	25.4 - 27.8 %
Nitrogen	9.9 - 10.8 %

Specification: PRD.1.ZQ5.10000022905

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.