


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
3-Hydroxy-4-nitrobenzaldehyde - 97%

Product Number: 156167
CAS Number: 704-13-2
MDL: MFCD00007109
Formula: C₇H₅NO₄
Formula Weight: 167.12 g/mol



| TEST | Specification |
|--|--------------------------|
| Appearance (Color) | Yellow to Brown |
| Appearance (Form) | Conforms to Requirements |
| Powder, Crystals, Crystalline Powder and/or Chunks | |
| Infrared spectrum | Conforms to Structure |
| Carbon | 48.5 - 52.1 % |
| Nitrogen | 8.0 - 8.7 % |
| Purity (GC) | ≥ 96.5 % |

Specification: PRD.0.ZQ5.10000038365

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.