

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
β-Cyclodextrin polymer - powderProduct Number: **C2485**

MDL: MFCD03456204

Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Light Yellow
Solubility (Turbidity)	Clear
100 mg plus 1 mL, H ₂ O	
Proton NMR spectra	Conforms to Structure
Solubility	Conforms
Solubility in water is > = 0.25 g/ml	

Specification: PRD.0.ZQ5.10000020682

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.