

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Magnesium bromide ethyl etherate - 99%

Product Number:

225959

CAS Number:

29858-07-9

MDL:

MFCD00064500

Formula:

Br₂Mg · C₄H₁₀O

Formula Weight:

258.23 g/mol



TEST

Specification

Appearance (Color)

White to Beige

Appearance (Form)

Powder or Solid or Chunk(s)

Infrared spectrum

Conforms to Structure

Titration by AgNO₃

98.5 - 101.5 %

Specification: PRD.2.ZQ5.10000013591

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.