

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Chloro(tri-tert-butylphosphine)gold(I)

Product Number:

**682705**

CAS Number:

69550-28-3

MDL:


MFCD09753040

Formula:

C<sub>12</sub>H<sub>27</sub>AuClP

Formula Weight:

434.74 g/mol



### TEST

### Specification

Appearance (Color)

White

Appearance (Form)

Powder or Crystals or Granules or Chunks

Proton NMR spectrum

Conforms to Structure

Phosphorus NMR Spectra

Conforms to Structure

Carbon

32.2 - 34.1 %

Specification: PRD.1.ZQ5.10000036624

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.