

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
mono-2-(Methacryloyloxy)ethyl maleate

Product Number: 483729
CAS Number: 51978-15-5
MDL: MFCD00080563
Formula: C₁₀H₁₂O₆
Formula Weight: 228.20 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	Colorless to Yellow
Appearance (Form)	Viscous Liquid
Might Contain Precipitate	
Proton NMR spectra	Conforms to Structure
Purity (Titration by NaOH)	≥ 80.0 %

Specification: PRD.1.ZQ5.10000014616

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.