


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Copper(II) acetate hydrate - 98%

Product Number: 341746
CAS Number: 66923-66-8
MDL: MFCD00149570
Formula: C₄H₆CuO₄ · xH₂O
Formula Weight: 181.63 g/mol



TEST	Specification
Appearance (Color)	Conforms to Requirements
Blue-Green to Dark Green	
Appearance (Form)	Crystals
X-Ray Diffraction	Conforms to Structure
Compleximetric EDTA	97.5 - 102.5 %
Water Corrected	
Moisture Content (LOD)	≥ 4.7 %

Specification: PRD.0.ZQ5.10000070315

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.