


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
3,5-Dimethyl-4-methoxyphenylboronic acid - ≥95%

Product Number: 598062
CAS Number: 301699-39-8
MDL: MFCD01114648
Formula: C₉H₁₃BO₃
Formula Weight: 180.01 g/mol



TEST

Specification

Appearance (Color)
Appearance (Form)
Solid, Powder, Crystals, Crystalline Powder, and/or Chunks
Proton NMR spectrum
Titration with NaOH

White to Yellow
Conforms to Requirements

Conforms to Structure
≥ 95.0 %

Specification: PRD.0.ZQ5.10000038015

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.