

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Calcium - dendritic pieces, purified by distillation, 99.5% trace metals basis

Product Number:

465895

CAS Number:

7440-70-2

MDL:

MFCD00085314

Formula:

Ca

Formula Weight:

40.08 g/mol

Ca

TEST	Specification
Appearance (Color)	Silver or Grey
Appearance (Form)	Chunk(s)
ICP Major Analysis	Confirmed
Confirms Calcium Component	
Purity	Meets Requirements
99.5% Based On Trace Metals Analysis	
Trace Metal Analysis	≤ 5500.0 ppm

Specification: PRD.0.ZQ5.10000043008

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.