

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com


Product Specification

Product Name:
Tannic acidProduct Number: **W304204**

CAS Number: 1401-55-4

Formula: C₇₆H₅₂O₄₆

Formula Weight: 1,701.2 g/mol

**TEST****Specification**

Appearance (Color)	Light Beige to Dark Beige to Dark Brown
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Carbon	45.5 - 61.6 %
Water (by Karl Fischer)	≤ 10.0 %
Arsenic (As)	≤ 3 ppm
Cadmium (Cd)	≤ 1 ppm
Mercury (Hg)	≤ 1 ppm
Lead (Pb)	≤ 10 ppm
Expiration Date Period	-----
60 Months	

Specification: PRD.3.ZQ5.10000076725

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.