

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com **Email USA:** techsery@sial.com Outside USA: eurtechserv@sial.com

## **Product Specification**

**Product Name:** 

Titanium - rod, diam. 6.35 mm, 99.99% trace metals basis

**Product Number:** 347132 Τi CAS Number: 7440-32-6

MFCD00011264 MDL:

Formula:

47.87 g/mol Formula Weight:

**TEST Specification** 

Appearance (Color) Conforms

Silver

Appearance (Form) Conforms

Rods

6.35 mm

Diameter Conforms

Confirmed ICP Major Analysis

Confirms Titanium Component

Trace Metal Analysis < 200.0 ppm

Purity Meets Requirements

99.99% Based On Trace Metals Analysis

Specification: PRD.0.ZQ5.10000026231