

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Nerve Growth Factor-2.5S from murine submaxillary gland - NGF-2.5S, lyophilized powder, suitable for cell culture

Product Number:

N6009

CAS Number:

93928-24-6

MDL:

MFCD00161772

Storage Temperature:

-20 °C

TEST	Specification
Appearance (Form)	Powder
Protein	Conforms
10 µg/vial (for 10 µg SKU)	
Purity by SDS Page	Conforms
> 95%	
Bioburden	
Report result	
Mycoplasma Test	None Detected

Specification: PRD.1.ZQ5.10000102902

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

