

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Silver nitrate – ACS reagent, ≥99.0%

Product Number: 209139  
CAS Number: 7761-88-8  
MDL: MFCD00003414  
Formula: AgNO<sub>3</sub>  
Formula Weight: 169.87 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Crystals
X-Ray Diffraction	Conforms to Structure
Titration by KSCN (% Purity)	≥ 99.0 %
Appearance (Clarity)	Pass
Chloride (CL)	≤ 5 ppm
Copper (Cu)	≤ 2 ppm
Iron (Fe)	≤ 2 ppm
Lead (Pb)	≤ 0.001 %
Sulfate	≤ 0.002 %
Free Acid	Pass
Substances not ppt. by HCl	≤ 0.01 %
Recommended Retest Period	-----
2 Years	
Meets ACS Requirements	Current ACS Specification

Specification: PRD.1.ZQ5.10000024086

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.