

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Vanadium(III) chloride - 97%

Product Number: 208272
CAS Number: 7718-98-1
MDL: MFCD00011454
Formula: Cl₃V
Formula Weight: 157.30 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Purple to Very Dark Purple and Dark Blue to Very Dark Blue and Black
Appearance (Form)	Powder
Titration with AgNO ₃	96.5 - 103.5 %
ICP Major Analysis	Confirmed
Confirms Vanadium Component	

Specification: PRD.3.ZQ5.10000092967

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.