

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com


Product Specification

Product Name:

Methanol solution - contains 0.10 % (v/v) formic acid, UHPLC, for mass spectrometry, ≥99.5%

Product Number:

632546



MDL:

MFCD00004595

Formula:

CH₄O

Formula Weight:

32.04 g/mol

TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
HPLC Absorbance 254 nm	≤ 10 mAU
HPLC Absorbance 210 nm	≤ 20 mAU
Purity (GC)	≥ 99.50 %A
Suitability	Suitable
LC-MS application tested	
UV Transmittance 230nm	≥ 10 %
UV Transmittance 254nm	≥ 90 %
Calcium (ICP)	≤ 0.5 ppm
Potassium (ICP)	≤ 0.5 ppm
Magnesium (ICP)	≤ 0.5 ppm
Sodium (ICP)	≤ 2 ppm
Expiration Date Period	-----
18 months	

Specification: PRD.5.ZQ5.10000017638

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.