

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Polyethylene - High density, melt index 12 g/10 min (190 °C/2.16kg)

Product Number:

427985


CAS Number:

9002-88-4

MDL:

MFCD00084423

Formula:

C₂H₄**TEST****Specification**

Appearance (Color)

White to Off-White

Appearance (Form)

Beads or Pellets

Infrared Spectrum

Conforms to Structure

Melt Index

≥ 10.5

g/10min (2.16 kg at 190°C)

Specification: PRD.1.ZQ5.10000038740

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.