


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Ethylene glycol monopropyl ether – ReagentPlus®, 99.4%

Product Number: 537675
CAS Number: 2807-30-9
MDL: MFCD00072645
Formula: C₅H₁₂O₂
Formula Weight: 104.15 g/mol



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Water (by Karl Fischer)	≤ 0.100 %
Specific Gravity At 20 Degrees Celsius	0.912 - 0.914
Acidity as Acetic Acid	≤ 0.01 %
Miscellaneous Supplier Data 2-Propoxyethanol	≥ 99.40 %
Color Test APHA	≤ 10

Specification: PRD.1.ZQ5.10000025330

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.