


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2,3-Dichloromaleic anhydride - 97%

Product Number: **D65003**
CAS Number: 1122-17-4
MDL: MFCD00005520
Formula: C₄Cl₂O₃
Formula Weight: 166.95 g/mol



TEST

Specification

Appearance (Color)

Conforms to Requirements

White to Beige to Light Yellow

Appearance (Form)

Conforms to Requirements

Crystalline Flakes, Powder or Chunks

Infrared spectrum

Conforms to Structure

Titration by AgNO₃

96.5 - 103.5 %

After Oxygen Combustion

Purity (GC)

≥ 96.5 %

Specification: PRD.0.ZQ5.10000026871

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.