

3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Dulbecco's Modified Eagle's Medium/Ham's Nutrient Mixture F12 – cell culture tested

Product Number: 51448C

MDL: MFCD00217417

Storage Temperature: 2 - 8 °C

TEST	Specification
Appearance (Turbidity)	Clear
Appearance (Color)	Red
Appearance (Form)	Solution
pH	7.0 - 7.4
Osmolality	290 - 330 mOs/kg
Sterility by USP guidelines	No Microbial Growth Detected
Endotoxin Level	≤ 1.0 EU/ml

Specification: PRD.0.ZQ5.10000037522

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.