


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Tris(2,2,6,6-tetramethyl-3,5-heptanedionato)manganese(III) - 97%

Product Number: **700223**
CAS Number: 14324-99-3
Formula: C₃₃H₅₇MnO₆
Formula Weight: 604.74 g/mol



TEST	Specification
Appearance (Color)	Very Dark Brown to Black
Appearance (Form)	Conforms to Requirements
Powder or Crystals	
Infrared spectrum	Conforms to Structure
Carbon	63.3 - 67.8 %
ICP Major Analysis	Confirmed
Confirms Manganese Component	

Specification: PRD.0.ZQ5.10000032093

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.