

3050 Spruce Street, Saint Louis, MO 63103, USA  
Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)  
Email USA: [techserv@sial.com](mailto:techserv@sial.com)  
Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

Product Specification

Product Name:  
Hydroquinone tablets

Product Number: HT1014  
CAS Number: 123-31-9

Formula Weight: 110.11 g/mol  
Storage Temperature: 2 - 8 °C

TEST	Specification
Appearance (Color)	White to Off White
Appearance (Form)	Round Tablet
Intact	
Performance	Conforms
Using Sigma-Aldrich Procedure HT101,spirochetes, helicobacter and other non-filamentous bacteria stain dark brown to black, background should be bright yellow to light brown	
Expiration Date Period	-----
24 months	

Specification: PRD.1.ZQ5.10000028016

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.