

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Copper(I) 3-methylsalicylate

Product Number: 715174

Formula: C₈H₇CuO₃
Formula Weight: 214.69 g/mol



TEST	Specification
------	---------------

Appearance (Color)	Off White to Light Green and Faint Brown to Brown and Faint Red to Red
Appearance (Form)	Powder or Solid
Infrared Spectrum	Conforms to Structure
Complexometric EDTA % Copper	28.6 - 30.6 %

Specification: PRD.2.ZQ5.10000023570

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

