

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

GlycoProfile™ IV Chemical Deglycosylation Kit - TFMS Deglycosylation System

Product Number:

PP0510

TEST	Specification
------	---------------

Suitability

Suitable

Not less than 95% deglycosylation of RNase B observed by SDS and judged by band shift from 15 to 13.5 kDa.

Specification: PRD.1.ZQ5.10000009298.000

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.