

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Monoclonal Anti-Netrin-1 antibody produced in mouse - clone 158936, purified immunoglobulin, lyophilized powder

Product Number: **N3787**

MDL: MFCD02099705

Storage Temperature: -20 °C

TEST	Specification
Formulation	Conforms
Lyophilized from 0.2 um-filtered solution in phosphate buffered saline containing carbohydrates	
Working Dilution	Conforms
by Direct Immunoblotting	
1-2 ug/mL detects mouse Netrin-1 at approx 25 ng/lane under non-reducing and reducing conditions	

Specification: PRD.2.ZQ5.10000070725

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.