

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-Synapsin I antibody produced in rabbit - affinity isolated antibody, lyophilized powder

Product Number: **S193**

Storage Temperature: -20 °C

TEST	Specification
Formulation	Conforms
Lyophilized from 5 mM ammonium bicarbonate	
Reconstitution	Conforms
Adequate for 10-mini-Western blots when reconstituted with 50 microliters PBS.	
Working Dilution	Conforms
by Western Blot	
1:1,000	
Working Dilution	Conforms
by Immunoprecipitation	
1 ug will quantitatively immunoprecipitate all synapsin in 200 ug rat lysate	

Specification: PRD.1.ZQ5.10000043666

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.