

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Potassium [(4-methoxybenzyloxy)methyl]trifluoroborate - 97%

Product Number:

767441

CAS Number:


1027642-26-7

Formula:

C₉H₁₁BF₃KO₂

Formula Weight:

258.09 g/mol



TEST

Specification

Appearance (Color)

White to Tan

Appearance (Form)

Powder or Chunks

Fluorine NMR Spectra

Conforms to Structure

Boron HNMR Spectra

Conforms to Structure

Proton NMR spectra

Conforms to Structure

Carbon

40.4 - 43.3 %

Purity Titration NAOH by HCL

96.5 - 103.5 %

Specification: PRD.0.ZQ5.10000027865

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.