

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Anti-phospho-GATA4 (pSer¹⁰⁵) antibody produced in rabbit - affinity isolated antibody, buffered aqueous solutionProduct Number: **G8044**

MDL: MFCD06798487

Storage Temperature: -70 °C

TEST	Specification
------	---------------

Formulation	Conforms
-------------	----------

Solution in Dulbecco's PBS, pH 7.3 with
1 mg/mL BSA and 0.05% Sodium Azide

Note	Conforms
------	----------

sufficient for 10 blots

Working Dilution	Conforms
------------------	----------

By Western Blot 1:500

Specification: PRD.1.ZQ5.1000088460

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.