

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Notch-2/Fc Chimera from rat – recombinant, expressed in NSO cells, ≥90% (SDS-PAGE), lyophilized powder

**Product Number:** N9161

MDL: MFCD06412282

Storage Temperature: -20 °C

TEST	Specification
Appearance (Form)	Lyophilized Powder
Formulation	Conforms
Lyophilized from 0.2 µm-filtered solution of phosphate buffered saline	
Titer by Functional Elisa	Conforms
5 µg/ml (100µl/well) can bind recombinant rat jagged-1/FC with a linear range of 6-400 ng/ml	
Endotoxin Assay	Accept
< 1.0 EU/µg cytokine	
Purity of SDS PAGE	Conforms
> 90%	

Specification: PRD.1.ZQ5.10000026596

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.