


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Benzophenone - ReagentPlus[®], 99%

Product Number: B9300
CAS Number: 119-61-9
MDL: MFCD00003076
Formula: C₁₃H₁₀O
Formula Weight: 182.22 g/mol



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Flakes
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 98.5 %
Solubility (Turbidity) 100mg/ml, Ethanol	Clear
Solubility (Color)	Colorless to Faint Yellow

Specification: PRD.0.ZQ5.10000043981

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.