


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
L-Ascorbic acid - plant cell culture tested

Product Number: A2174
CAS Number: 50-81-7
MDL: MFCD00064328
Formula: C₆H₈O₆
Formula Weight: 176.12 g/mol



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
50 mg/mL, H ₂ O	
Titration with Iodine	≥ 98 %
Plant Cell Culture Test	Pass

Specification: PRD.0.ZQ5.10000028465

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.