

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Glass beads, acid-washed - 710-1,180 µm (16-25 U.S. sieve)

Product Number:

G1152

TEST	Specification
------	---------------

Appearance (Color)

Colorless to Faint Green to Faint Yellow

Appearance (Form)

Beads

Appearance (Clarity)

Clear

Particle Size Distribution

≥ 90 %

Within 710 to 1180 microns

Recommended Retest Period

5 years

Specification: PRD.0.ZQ5.10000045142

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.