

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Formaldehyde-sodium bisulfite adduct - 95%

Product Number:

112704

CAS Number:

870-72-4

MDL:


MFCD00044664

Formula:

CH<sub>3</sub>NaO<sub>4</sub>S

Formula Weight:

134.09 g/mol

**TEST****Specification**

Appearance (Color)

White

Appearance (Form)

Powder

Infrared spectrum

Conforms to Structure

Titration with Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub>

94.5 - 105.5 %

Specification: PRD.1.ZQ5.10000029479

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.