


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Water – PCR Reagent

Product Number: **W1754**  
CAS Number: 7732-18-5  
MDL: MFCD00011332  
Formula: H<sub>2</sub>O  
Formula Weight: 18.02 g/mol



TEST	Specification
------	---------------

Appearance (Turbidity)

Clear

Appearance (Color)

Colorless

Appearance (Form)

Liquid

DNase, Exonuclease Detection

None Detected

NICKase, Endonuclease Detection

None Detected

RNAse Detection

None Detected

Suitable for PCR

Suitable

Specification: PRD.0.ZQ5.10000063815

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.