

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Potassium 2-methoxy-3-(trifluoromethyl)pyridine-4-trifluoroborate - 95%

Product Number: 777382

Formula: C₇H₅BF₆KNO
Formula Weight: 283.02 g/mol



TEST

Specifi

Appearance (Color)	White to Light Yellow and Faint Brown to Light Brown
Appearance (Form)	Powder or Crystal or Chunk(s)
Fluorine NMR Spectra	Conforms to Structure
Boron HNMR Spectra	Conforms to Structure
Proton NMR Spectrum	Conforms to Structure
Carbon	28.1 - 31.3 %
Nitrogen	4.7 - 5.2 %
Titration with HClO ₄	94.5 - 105.5 %

Specification: PRD.0.ZQ5.10000039285

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.