

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Hyperbranched bis-MPA polyester-32-hydroxyl, generation 3 - ≥97%

Product Number: **686581**
CAS Number: 326794-48-3

Formula: C₁₅₅H₂₅₆O₉₃
Formula Weight: 3,607.64 g/mol
Storage Temperature: 2 - 8 °C

TEST	Specification
Appearance (Color)	White
Appearance (Form) Powder and/or Chunks	Conforms to Requirements
Infrared spectrum	Conforms to Structure
Carbon	49.8 - 52.9 %
Molecular Number by GPC Specification: Report Results	
Molecular Weight by GPC Specification: Report Results	
Poly Dispersity MW/MN Specification: Report Results	
Specification: PRD.0.ZQ5.10000027349	

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.