

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Potassium hexachloroplatinate(IV) - ≥99.9% trace metals basis

Product Number: 520861
CAS Number: 16921-30-5
MDL: MFCD00011389
Formula: Cl₆K₂Pt
Formula Weight: 485.99 g/mol



TEST	Specification
Appearance (Color)	Yellow to Yellow-Orange
Appearance (Form)	Powder
ICP Major Analysis Confirms K and Pt Components	Confirmed
Gravimetric Analysis % Pt	38 - 41 %
Trace Metal Analysis	≤ 1000 ppm
Purity ≥99.9% Based on Trace Metals Analysis	Meets Requirements

Specification: PRD.1.ZQ5.1000008887

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.