


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
4-Chloro-1-naphthol - tablet

**Product Number:** C6788  
**CAS Number:** 604-44-4  
**MDL:** MFCD00003974  
**Formula:** C<sub>10</sub>H<sub>7</sub>ClO  
**Formula Weight:** 178.61 g/mol  
**Storage Temperature:** -20 °C



TEST	Specification
Appearance	Conforms
Off-White to Tan	
Diameter	5/32 inches
Form	Round Tablet
Dissolution Time	≤ 3 min
Solubility (Color)	Colorless
Solubility (Turbidity)	Clear
Solubility (Form)	Solution
Total Tablet Weight	28.8 - 35.2 mg/tab
Dot Blot	Conforms
Compares to untableted substrate	
Note	-----
Dissolution time and solubility are the results of one tablet dissolved in 10 mL of MeOH	

Specification: PRD.1.ZQ5.10000008046

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.