


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Flurbiprofen - cyclooxygenase inhibitor

Product Number: F8514
CAS Number: 5104-49-4
MDL: MFCD00079303
Formula: C₁₅H₁₃FO₂
Formula Weight: 244.26 g/mol



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity)	Clear
50 mg/ml, MeOH	
Carbon	72.3 - 74.8 %
Purity (HPLC)	≥ 98.5 %
Recommended Retest Period	-----
3 years	

Specification: PRD.1.ZQ5.10000026242

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.