


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Glycine sodium salt hydrate - 98%

Product Number: 219517
CAS Number: 207300-76-3
MDL: MFCD00150683
Formula: C₂H₄NO₂Na · xH₂O
Formula Weight: 97.05 g/mol



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder or Chunks
Infrared spectrum	Conforms to Structure
Titration with H ₂ SO ₄ Dried Basis	97.5 - 102.5 %
Water (by Karl Fischer)	≤ 24.0 %

Specification: PRD.1.ZQ5.10000024290

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.