


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2-Ketohexanoic acid sodium salt - 97.0-103.0%

Product Number: K6625
CAS Number: 13022-85-0
MDL: MFCD00013063
Formula: C₆H₉NaO₃
Formula Weight: 152.12 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Solubility (Color)	Colorless
Solubility (Turbidity) At 50 mg/ml of H ₂ O	Clear
Proton NMR spectrum	Conforms to Structure
Titration with HClO ₄	97.0 - 103.0 %
Purity (TLC)	≥ 97 %

Specification: PRD.2.ZQ5.10000009122

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.