


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Ketorolac tris salt - ≥99%, crystalline

Product Number: K1136
CAS Number: 74103-07-4
MDL: MFCD00887595
Formula: C₁₅H₁₃NO₃ · C₄H₁₁NO₃
Formula Weight: 376.40 g/mol



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Powder
HNMR	Conforms to Structure
Solubility (Color)	Colorless to Faint Yellow
Solubility (Turbidity)	Clear
15 mg/mL, H ₂ O	
Carbon	59.7 - 61.5 %
Nitrogen	7.1 - 7.7 %
Purity (TLC)	≥ 99 %

Specification: PRD.0.ZQ5.10000098273

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.