

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Potassium fluoride – anhydrous, powder, ≥99.9% trace metals basis

Product Number: 449148
CAS Number: 7789-23-3
MDL: MFCD00011398
Formula: FK
Formula Weight: 58.10 g/mol

KF

TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
ICP Major Analysis Confirms K Component	Confirmed
Purity > = 99.9% Based On Trace Metals Analysis	Conforms
Trace Metal Analysis	≤ 1000.0 ppm

Specification: PRD.0.ZQ5.10000059892

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.