


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
9-Anthracenecarboxaldehyde - 97%

Product Number: 278688
CAS Number: 642-31-9
MDL: MFCD00001254
Formula: C₁₅H₁₀O
Formula Weight: 206.24 g/mol



TEST

Specification

Appearance (Color)

Yellow to Dark Yellow and Yellow-Brown
and Yellow-Green

Appearance (Form)

Conforms to Requirements

Powder or Crystalline Powder

Infrared spectrum

Conforms to Structure

Purity (HPLC)

≥ 96.5 %

Specification Date : 04/15/2011

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.