


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
3-Amino-2-methylpyridine - 97%

Product Number: 662690
CAS Number: 3430-10-2
MDL: MFCD03788195
Formula: C₆H₈N₂
Formula Weight: 108.14 g/mol



TEST

Specification

Appearance (Color)

White to Yellow and Faint Beige to Beige
and Faint Brown to Light Brown

Appearance (Form)

Conforms to Requirements

Powder, Crystals or Crystalline Powder and/or Chunks

Infrared spectrum

Conforms to Structure

Purity (Titration by HCL04)

96.5 - 103.5 %

Purity (GC)

≥ 96.5 %

Specification: PRD.1.ZQ5.10000005824

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.